



LV Rescue Kit

kit Instructions:

- Before work commences, the presence and correct assembly of all components shall be verified and each component suitably inspected for signs of damage, deterioration and legibility of markings.
- Rescue kits and their components shall not be altered in any way that may affect their operation as originally recommended or supplied.

Periodic Inspections:

- Every 6 months The Rescue Kit must be examined. Any component failing test or showing signs of wear or deterioration, that would effect its safe operation must be replaced.



Overview / Features

Item No	Name	Description	Qty	Image
26008	Kit Bag	460mm x 120mm x 240mm, Orange PVC Bag with Velcro Fastening + CPR Card		
GAN-CRK	LV Crook	370mm Long x 30mm Dia, Grey PVC, Tested 5KV	1 EA	
CAT6L2D	LED Torch	Camelion Yellow Torch, 2 x D Batts, Non Conductive	1 EA	
BAR-ISO	"Isolate Here" Tag	260 x 180 x 1.5mm, White HDPE With Text, CNV Rope	1 EA	
WFB1812	Fire Blanket	1.8m Long x 1.2m Wide	1 EA	
TYC-MTD	Multi Trauma Dressing	30" Long x 10" Wide	1 EA	
AE0650360	LV Gloves	650 volt 360mm Long Gloves	1 Pr	



WEB: www.austutilities.com.au
 EMAIL: sales@austutilities.com.au

PHONE: 08 9406 3200